

ElectroEB10B Competent Cells

Product	Cat. No.	Transformations
ElectroEB10B Competent Cells	12573	30

Description

ElectroEB10B Competent Cells have been manufactured using a proprietary technology rendering the cells highly efficient for DNA uptake when used in an electroporation device. To utilize the cells at their highest efficiency, a recommended transformation protocol is included with each kit.

ElectroEB10B Competent Cells are only provided in tubes. All kits include a test plasmid for quality control purposes. Each tube contains pre-dispensed cells in 100 μ l aliquots, enough for 5 transformations using the recommended transformation protocol. Full processing time, including recovery, is about 1 hour to ensure the highest level of transformation. Edge BioSystems guarantees transformation efficiencies exceeding 10¹⁰ cfu/ μ g pUC19.

Genotype: F⁻ mcrA Δ (*mmr*-*hsdRMS*-*mcrBC*) Φ 80*lacZ* Δ M15 Δ *lacX74 recA1 endA1 araD139 Δ (*ara, leu*)7697 *galU galK* λ ⁻ *rpsL nupG tonA*.*

Kit Components	12573
ElectroEB10B Competent Cells	6 tubes
pUC19 Supercoiled DNA, 100ng/ml	1 tube

Quality Control

Each lot has been tested to assure high transformation efficiency using 10pg pUC19 supercoiled DNA and the recommended protocol. Transformation efficiency will exceed >10¹⁰ cfu/ μ g pUC19 under these conditions.

Equipment and Materials Not Provided

- SOC medium for recovery: 20g/l tryptone, 5g/l yeast extract, 10mM NaCl, 2.5mM KCl, 10mM MgCl₂, 10mM MgSO₄, 20mM glucose (MgCl₂, MgSO₄ and glucose should be added after autoclaving)
- An orbital shaker
- An electroporation device with 1mm cuvettes
- 14ml round bottom culture tubes (1 tube per transformation)

- LB-agar selective plates or selective liquid medium
- Pipettors

Storage Conditions

ElectroEB10B Competent Cells should be stored in a -80°C freezer. Please note that competent cells are very sensitive to cycles of freezing and thawing and should not be exposed to temperature variations.

Recommended Protocol for Tubes

- Place tubes in ice and wait approximately 5-10 minutes until they thaw.
- Chill 1 microcentrifuge tube and 1 electroporation cuvette on ice for each sample.
- When thawed, mix cells by gently tapping the bottom of the tube. Aliquot 20-40 μ l of cells into each chilled microcentrifuge tube. Control reaction: Use 20 μ l of cells for the control reaction. Excess cells can be re-frozen in a dry ice/ethanol bath and returned to the -80°C freezer.
- Pipette the DNA to be transformed (1-5 μ l) to the bottom of the tube and mix by gently tapping the bottom of the tube. Control transformation: Dilute pUC19 supercoiled DNA 1:10 with sterile H₂O, then add 1 μ l of the diluted pUC19 supercoiled DNA to one of the tubes. Discard diluted pUC19 supercoiled DNA after use.

Note: Refrozen cells will have reduced transformation efficiency.

Note: Do not mix by pipetting up and down since that will lower the transformation efficiency. Mix and electroporate only 1 sample at a time.

- Immediately transfer cell/DNA mix into a pre-chilled electroporation cuvette.
- Electroporate samples using the following conditions for the BioRad GenePulser® II: 2.0kV, 200 Ω , 25 μ F.
- Add 1 ml SOC medium to the cuvette and transfer the sample to a 14ml round bottom tube.
- Transfer the tubes to a shaking incubator and then shake at 300 rpm at 37°C for 1 hour.

Note: Do not introduce bubbles into the sample.

9. Dilute samples as necessary with SOC medium. Control reaction: Dilute the pUC19 control reaction 1:100 with S.O.C. medium.
10. Spread 100-200µl of sample on pre-warmed LB-agar selective plates or inoculate into selective liquid medium. Control transformation: Spread 50µl on LB-ampicillin agar plates and expect 50 colonies ($>10^{10}$ cfu/µg pUC19).

Notes:

- a. Transformation efficiencies for ligation mixtures will be 10-100 fold lower than pUC19 supercoiled DNA ligations. It is important to remove salts and buffers from any samples prior to electrotransformation since they will inhibit electroporation. Ligation reactions should be ethanol precipitated or diluted before electrotransformation. Ligations can be diluted 1:5 with dH₂O. Excess salt, buffer or DNA can cause electric arcing.
- b. Using an electroporator other than the BioRad GenePulser® II may require setting adjustments in order to obtain the optimal efficiency for that specific model.
- c. Calculating transformation efficiency

$$\frac{(\text{cfu on control plate})}{(\text{pg of supercoiled vector})} \times (10^6 \text{ pg} / \mu\text{g}) \times (\text{final dilution}) = \text{cfu} / \mu\text{g DNA}$$

(Note: cfu = colony forming units)

For example, if the control reaction plate has 50 colonies when 50µl of diluted reaction is plated then,

$$\frac{(50 \text{ cfu})}{(10 \text{ pg})} \times \frac{(10^6 \text{ pg})}{(\mu\text{g})} \times (1 \text{ ml}) / (0.05 \text{ ml plated}) \times 10^2 = 1.0 \times 10^{10} \text{ cfu}/\mu\text{g pUC19 DNA}$$