

EB5Alpha Competent Cells

Product	Cat. No.	Transformations
EB5Alpha Competent Cells, Tubes	83254	12

Description

EB5Alpha Competent Cells are chemically competent and have been manufactured using a proprietary technology that renders the cells highly efficient for DNA uptake. To utilize the cells at their highest efficiency, a recommended transformation protocol is included with each kit.

EB5Alpha Competent Cells are available in single use tubes that provide a simple and reliable method for high-efficiency, single use transformation. All kits include a test plasmid for quality control purposes. Cells are pre-dispensed in 50 μ l aliquots (tubes).

Full processing time for 1 tube (including recovery) is about 1 hour 20 minutes to ensure the highest level of transformation. Edge BioSystems guarantees transformation efficiencies that exceed 10⁹ cfu/ μ g pUC19.

Genotype: F⁻ Φ 80*lacZ* Δ M15 Δ (*lacZYA-argF*) U169 *recA1 endA1 hsdR17*(r_k⁻,m_k⁺) *phoA supE44 thi-1 gyrA relA1 tonA*.

- 14ml round bottom culture tubes (1 tube per single use aliquot).
- LB-agar selective plates or selective liquid medium.

Storage Conditions

EB5Alpha Competent Cells should be stored in a -80°C freezer. Please note that competent cells are very sensitive to cycles of freezing and thawing and should not be exposed to temperature variations.

Recommended Protocol for Tubes

- Place single use tubes on ice and wait approximately 5 minutes until they thaw.
- Pipette the DNA to be transformed to the bottom of the tube and mix by pipetting 50 μ l of air to the bottom of the tube. Control transformation: Dilute pUC19 supercoiled DNA 1:10 with sterile H₂O, then add 1 μ l of the diluted pUC19 supercoiled DNA to one of the tubes. Discard diluted pUC19 supercoiled DNA after use.

Note: Do not mix by pipetting up and down since that will lower the transformation efficiency.

- Incubate the tubes in ice for 10 minutes.
- Transfer the tubes to a 42°C water bath, incubate for 40 seconds and transfer the tubes back to ice.
- Incubate the tubes for 2 minutes in ice.
- Transfer the cells into a 14ml round bottom culture tubes filled with 1ml of pre-warmed SOC medium and then shake at 300 rpm at 37°C for 1 hour.
- Plate on pre-warmed LB-agar selective plates or inoculate into selective liquid medium. For the control transformation with pUC19 supercoiled DNA, plate 10 μ l on LB-ampicillin agar plates and expect 100 colonies (>10⁹ cfu/ μ g pUC19).

Kit Components 83254

EB5Alpha Competent Cells	12 tubes
pUC19 Supercoiled DNA, 100ng/ml	1 tube

Quality Control

Each lot has been tested to assure high transformation efficiency using 10pg pUC19 supercoiled DNA and the recommended protocol. Transformation efficiency will exceed >10⁹ cfu/ μ g pUC19 under these conditions.

Equipment and Materials Not Provided

- SOC medium for recovery: 20g/l tryptone, 5g/l yeast extract, 10mM NaCl, 2.5mM KCl, 10mM MgCl₂, 10mM MgSO₄, 20mM glucose (MgCl₂, MgSO₄ and glucose should be added after autoclaving).
- An orbital shaker.
- A 42°C water bath.

Warning: This product is intended for **research use only**. It is not to be used for diagnostic purposes in humans or animals.