



Troubleshooting Guide for MagDTR™ Technology

PROBLEM	POSSIBLE CAUSE	SUGGESTED SOLUTIONS
	<p>Incomplete aspiration of ethanol after binding and/or wash steps</p>	<ul style="list-style-type: none"> • Complete removal of ethanol is critical for successful cleanup. Residual ethanol can be removed by aspiration at a pipet setting of 150µl or higher, coupled with slow aspiration. There should be no need for a second aspiration. However, if a second aspiration is necessary, care should be taken to ensure that ethanol is not reintroduced into the wells from the pipet tip. • Increase the number of washes from 2 to 3 washes in order to ensure that all the salt and dye have been removed.
<p>Dye peaks present</p>	<p>Excessive dye present</p>	<ul style="list-style-type: none"> • Increase the number of washes from 2 to 3 washes to ensure that all salt and dye have been removed. • More than 4µl BigDye® used per sequencing reaction. • Use lower dye concentration. Dilute the dye with the appropriate sequencing buffer. • Increase the number of cycles during thermal cycling. • Increase the concentration of the template in the sequencing reaction.
	<p>Sequencing reaction conditions not optimal</p>	<ul style="list-style-type: none"> • Optimize template and primer concentrations. Consult dye manufacturer's recommendations. • Start with high purity template. • Check thermal cycler for program entry error. • Check for potential machine failure

Troubleshooting Guide for MagDTR™ Technology

PROBLEM	POSSIBLE CAUSE	SUGGESTED SOLUTIONS
Dye peaks present <i>cont...</i>	Wrong volume of ethanol used for precipitation	Use the correct volume of 100% ethanol for the reaction. Too high of a final concentration can lead to dye terminator retention.
	Wrong concentration of ethanol used for washes	<ul style="list-style-type: none"> • Check ethanol concentrations used. Be certain that freshly prepared 80% ethanol is used for the washes • Use at least 100µl of freshly prepared 80% ethanol per wash.
	Incompatible magnet	<ul style="list-style-type: none"> • Use the Edge BioSystems' MagWell Magnetic Separator. • Optimize the protocol for use with an alternative magnet.
	Incorrect storage of pre and post purified samples	<ul style="list-style-type: none"> • Avoid prolonged exposure to light. • Avoid extended storage at temperatures above 4°C. • Store post purified samples at -20°C.
Low recovery	Over dried samples	<ul style="list-style-type: none"> • Aspirate all ethanol. Reduce the drying time at room temperature to 3 minutes. Resuspend immediately • Aspirate all ethanol. Heat to 98°C for 30 seconds. Resuspend immediately. • At the elution step, add recommended volume of deionized water. Incubate for one minute at room temperature, then pipet mix to resuspend.

Troubleshooting Guide for MagDTR™ Technology

PROBLEM	POSSIBLE CAUSE	SUGGESTED SOLUTIONS
<p>Low recovery cont...</p>	<p>Incomplete elution at the resuspension step</p>	<ul style="list-style-type: none"> • Use the recommended volume of deionized water. Incubate for one minute at room temperature, then pipet mix to resuspend. • Vortex sample at medium speed for at least 30 seconds to properly resuspend the samples. • In the absence of a vortex, pipet aggressively.
	<p>Formamide used as an elution buffer</p>	<ul style="list-style-type: none"> • Direct elution with formamide is not recommended • Use the recommended volume of deionized water. Incubate for one minute at room temperature, then pipet mix to resuspend. • If the sample must be in formamide, the sample can be eluted in water then dried down and resuspended in formamide.
	<p>Incomplete drying of sample before resuspension</p>	<ul style="list-style-type: none"> • Increase drying time at room temperature or incubate at 98°C for 30 seconds.
	<p>Failure to precipitate sequencing products</p>	<ul style="list-style-type: none"> • Check ethanol concentration used. Be certain that 100% ethanol is used for precipitation. • Be certain that the correct volume of 100% ethanol was used.
	<p>Failure to retain sequencing products during washes</p>	<ul style="list-style-type: none"> • Check ethanol concentration used. Be certain to that freshly prepared 80% ethanol was used during the washes. • Try to remain within the suggested number of washes unless your reaction was done with greater than 4µl of BigDye®, then up to 3 washes may be required.

Troubleshooting Guide for MagDTR™ Technology

PROBLEM	POSSIBLE CAUSE	SUGGESTED SOLUTIONS
Low recovery cont...	Incomplete removal of reaction salts during the wash steps causing high background	<ul style="list-style-type: none"> Check ethanol concentration used. Be certain that freshly prepared 80% ethanol is used for the washes. Increases the number of washes from 2 to 3 washes to ensure that all salts and dye have been removed.
	Insufficient quantity of template used	Consult dye manufacturer's recommendations.
	Poor quality template used	<ul style="list-style-type: none"> Start with high purity template For optimal performance, use Edge BioSystems' QuickStep™2 for PCR purifications.
	Incompatible magnet used	<ul style="list-style-type: none"> Used Edge BioSystems' MagWell Magnetic Separator. Optimize protocol for use with an alternative method.
Noisy background on the sequence	Incomplete removal of reaction salts	<ul style="list-style-type: none"> Check ethanol concentration used. Be certain that freshly prepared 80% is used for the washes. Increase the number of washes from 2 to 3 washes to ensure that all salt and dye have been removed.
	Low signal strength or high signal strength	<ul style="list-style-type: none"> Follow the steps recommended for troubleshooting low signal sequences. If the signal is high, dilute the sample before loading onto the sequencer
	Incomplete drying of samples before resuspension	Increase drying time at room temperature or incubate at 98°C for 30 seconds.
	Contaminated template or primer	A clean template and a clean primer are required to obtain optimal results.



Troubleshooting Guide for MagDTR™ Technology

PROBLEM	POSSIBLE CAUSE	SUGGESTED SOLUTIONS
Noisy background on the sequence <i>cont...</i>	Thermal cycler failure	<ul style="list-style-type: none"> • Check thermal cycler program for entry error. • Check thermal cycler for mechanical or electrical failure. • Troubleshoot your thermal cycler to ensure optimal conditions for preparing reactions.
Resin settles too fast	Nature of the resin	<ul style="list-style-type: none"> • Shake well or vortex the resin prior to use. • If using a robotic workstation, make sure to pipet mix the resin prior to each aspiration step. • If using reagent reservoir, rock back and forth 8-10 times to resuspend resin before pipetting.
Resin forms aggregates during storage	Resin is exposed to strong magnetic field	Shake well or vortex the resin, then sonicate for 5-10 minutes before using.
	Resin is stored for extended period of time	Shake well or vortex the resin, then sonicate for 5-10 minutes before using.

For additional troubleshooting assistance, please contact Edge BioSystems at:

Telephone: (800) 326-2685 or (301) 990-2685	Fax: (301) 990-0881
Email: info@edgebio.com	Web: http://www.edgebio.com