



Troubleshooting Guide for PERFORMA[®] DTR 384-Well Plates

PROBLEM	POSSIBLE CAUSE	SUGGESTED SOLUTIONS
Signal intensity low	Centrifuge conditions were not optimal	Recalculate rpm to achieve selected speed for your centrifuge and reaction volume
	Sequencing reaction conditions were not optimal	Adjust sequencing reaction parameters; consult sequencing reagent manufacturer's recommendations
	Sample volume too low for sequencer	Add deionized water or sample buffer to purified reaction before loading on the sequencing instrument
Dye blobs	Too large a volume was loaded	<ul style="list-style-type: none"> • Load no more than 10 μl • If using > 2 μl BigDye[®] v 3.1 per 10 μl reaction, consider loading only 5 μl on the Performa[®] DTR 384-Well Plate • Consider using Performa[®] DTR Standard Plate for sample volumes > 10 μl, or sequencing reactions containing more than 3 μl BigDye[®]
	The concentration of BigDye [®] v. 3.1 was too high (<i>i.e.</i> More than 2 μ l BigDye [®] v 3.1 per 10 μ l reaction was loaded)	If the sequencing reaction conditions cannot be altered, consider loading only a portion of the reaction for cleanup and sequencing.
	Centrifuge conditions were not optimal	<ul style="list-style-type: none"> • Be certain that the pre-spin is for 2 min at 850 x g • Recalculate rpm to achieve selected speed for your centrifuge and reaction volumes
	Samples were not loaded properly	Load the samples slowly, drop-wise to the center of the gel matrix
	Signal intensity too high	Apply only a portion of sequencing reaction to the plate; high signal strengths exaggerate the effect of small quantities of dye terminator in the purified reaction



Troubleshooting Guide for PERFORMA[®] DTR 384-Well Plates

PROBLEM	POSSIBLE CAUSE	SUGGESTED SOLUTIONS
Dye blobs cont...	Sequencing reaction conditions were not optimal	Adjust sequencing reaction parameters; consult sequencing reagent manufacturer's recommendations
Truncation of 5' sequence	Centrifuge conditions were not optimal	Recalculate rpm to achieve selected speed for your centrifuge and reaction volumes
	Dye blobs present	See Dye blobs above
	Less than 5 μ l of sequencing reaction loaded	Adjust the sequencing reaction to 5 μ l before loading on the plate
	Sequencing reaction conditions not optimal	Adjust sequencing reaction parameters; consult sequencing reagent manufacturer's recommendations
Recovery volume higher than expected	Centrifuge conditions were not optimal	<ul style="list-style-type: none"> Be certain that the pre-spin is performed for 2 min at 850 x g Recalculate rpm for your centrifuge conditions
	More than 10 μ l of sequencing reaction loaded	Decrease volume of sequencing reaction loaded to 5-10 μ l
	Plate was stored incorrectly. If you observe a high recovery volume after the first spin, it is likely the plate has been frozen	<ul style="list-style-type: none"> Store sealed plates at 4°C Do not allow plates to freeze
Recovery volume lower than expected	Centrifuge conditions were not optimal	Recalculate rpm for your centrifuge and reaction volume
	Less than 5 μ l of sequencing reaction loaded	Increase volume of sequencing reaction to 5 -10 μ l
	Plates were stored incorrectly. If you observe a low recovery volume after the first spin, it is likely the plate has been dehydrated	<ul style="list-style-type: none"> Store plates at 4°C Plates sealed in original bag can be left at room temperature for short periods of time



Troubleshooting Guide for PERFORMA[®] DTR 384-Well Plates

For additional troubleshooting assistance, please contact Edge BioSystems at:

Telephone: (800) 326-2685 or (301) 990-2685	Fax: (301) 990-0881
Email: info@edgebio.com	Web: http://www.edgebio.com