


pEAK12 Vector Map



Sequences flanking the cloning site of pEAK12:

T7 promoter Hind III EcoR I
 5'-ATTAAT ACGACTCACTATAGGG AGACT TCTTTCTCCCATTTTCAGGTGTCGT AAGCTT GAATTC AA
Ase I

I-Ceu I VDE (PI-Sce I)
TAACTATAACGGTCCTAAGGTAGCGA AG {-5'- cDNA -3'-} TTTCATTACCTC TTTCTCCGCACCCACATA

AA CCTGCAGG CTCCCTTTAGTGAGGGTAATTCGTATGTGTATGATA GCGGCCGC-3'
Sse8387 I T3 promoter Not I